



Cyclic Voltammetry SOP

1. Set up should include the following:
 - a. Flask or vial with room for electrodes
 - b. Working electrode (Glassy carbon or other) attached to green wire
 - c. Counter Electrode (Pt wire) attached to red wire
 - d. Reference electrode (Ag) attached to the white wire
 - e. Electrolyte solution (2-5 ml, enough to cover electrodes) such as 0.1 M $[\text{NBu}_4][\text{PF}_6]$
2. Bubble inert gas such as Nitrogen into solution to clear Oxygen
3. First run the CV in the electrolyte to establish a background
4. Go to "T" for Technique and select CV
5. Adjust parameters:
 - Range will go from Init E (V) to final E(V)
 - Initial Scan polarity can be negative or positive depending on the direction of the Scan
 - Scan rate (V/s) is normally set at 0.1 (adjust accordingly for specific experiment)
 - Sweep segments 2 (every 2 is a cycle)
 - Sample interval (V) is normally 0.001 or can be adjusted accordingly
 - Quiet Time (sec) is normally 2
 - Sensitivity (A/V) is normally 1 e-004 or can be adjusted accordingly



6. Set Resistance Compensation.
 - a. By pressing "Test" in iR Compensation, Resistance should be <100.
 - b. If <100,
 - i. Move electrodes closer together
 - ii. Check electrode connections
 - iii. Add more electrolyte solution
 - c. iR compensation "Once" or Always
 - d. iR Comp Mode Automatic or manual and set resistance.



7. Run play button to take background
8. Add material (2-5 ml volume, 1.0-0.25 mmolar concentration, high concentration better)
9. With material, repeat step 6.
10. Run play button

11. Check Data Plot

- a. E_p Segment 1 = E_{p1}
- b. E_p Segment 2 = E_{p2}
- c. $E_{1/2} = (E_{p1} + E_{p2})/2$

12. For Electrosynthesis, Go to "T" for Technique and select electrosynthesis

13. In parameters set:

- a. $E_i = E_{1/2}$ from CV measurement
- b. Set Hold time till conversion